Optional: Add UCSC CpG island annotations

Partek Genomics Suite software can view annotation .BED files as tracks in the Genome Viewer. We can add a CpG islands track to the Genome Viewer using the UCSC Genome Browser CpG islands annotation.

- Go to UCSC Genome Browser
- Select Table Browser under Tools in the main command bar of the webpage (Figure 1)



Figure 1. Navigating to the Table Browser at the UCSC Genome Browser website

• Configure the Table Browser page as shown (Figure 2)

Â	Genomes	Genome Browser	Tools	Mirrors	Downloads	My Data	Projects	Help	About Us
Table B	rowser								
Use this For help queries. enrichm downloa clade: [group: table:] region: identified filter:] intersec correlat output fille type	program to re in using this a For more com mats, send the ded in their er Mammal Regulation pgIslandExt e genome (ers (names/ac create ction: create format: BED- file: cpg.bed e returned: (trieve the data asso application see <u>Using</u> applex queries, you m data to <u>GREAT</u> . Ref titrety from the <u>Sequ</u> genome: [Human 	ciated with g the Table iay want to fer to the 9 ience and track: C a gions p st upload a (lea compress	n a track in <u>a Browser fro</u> <u>b use Galax</u> <u>Credits pag</u> <u>Annotation</u> <u>y</u> <u>pG Islands</u> <u>pG Islands</u> <u>position chr</u> <u>list</u> <u>v</u> Sen ve blank to ed	text format, to ci or a description y or our <u>public 1</u> for the list of c <u>Downloads</u> pag assembly: Fe x :15,578,261-15,63 X:15,578,261-15,63 d output to keep output in b	alculate inters of the control <u>MySQL serve</u> ontributors an le. <u>ab. 2009 (GRCh</u> <u>d custom tracks</u> <u>21,068</u>	sections betw s in this form r. To examine ad usage rest 37/hg19) V track hubs lookup defir	veen track: , and the <u>i</u> e the biolog trictions as	s, and to retrieve DNA sequence covered by a track. <u>User's Guide</u> for general information and sample gical function of your set through annotation ssociated with these data. All tables can be
get outp	all user cart s	atistics	istom trac	ks) click he	re				
10 10301	un user call s	orango (including cu	Storn trac	(3), <u>click He</u>	<u>10</u> .				

Figure 2. Configuring the Table Browser to output CpG Islands BED file

- Set assembly to Feb. 2009 (GRCh37/hg19)
- Set group to Regulation
- Set *track* to CpG Islands
- Set table to cpgIslandExt
- Set *output format* to BED
- Set *output file* to **cpg.bed**
- Select get output

The Output cpglslandExt as BED page will open.

- Select get BED to download a compressed folder containing the BED file
- Unzip the file using 7-Zip, WinRAR, or a similar program of your choice to a location you will be able to find

Next, we can import the BED file into Partek Genomics Suite.

• Select Genomic Database... under Import under File in the main toolbar in Partek Genomics Suite (Figure 3)

🤣 Partek Genomics Suite - fen	nale_only ((female	_only)										- 🗆 X
File Edit Transform View	Stat Filt	ter To	ols Windo	w Custom Help									
Open													Workflows Methylation
Recent Files	I											i i	llumina BeadArray Methylation $\qquad imes$
Open As Child	Tout				0								Import
Import	Affyme	csv .txt, etrix) •		•								04/00
Open Project	Agilent		SM3453105 200492200025 D04C01							. K	Analusia		
Save Project Export Zipped Project	Two-co	olor Micro	oarray 🕨 🕨	GSM2452106_2004832000	2.	3.	4.	5.	6.	7.	8.	ΩK	Allalysis
Export Zipped Project	Illumina	a	•		Cell Type	Gender	cg09835024	cg14361672	cg12950382	cg02115394	cg12480	R	Visualization
Create New Spreadsheet	Genom	ic Datab	abase	_200483200025_R04C01	B cells	Female	-5.13333	3.88696	2.53662	-3.09657	-5.3429	2	Biological Interpretation
merge spreadsneets	Nimble	NimbleGen		_200483200021_R01C01	B cells	Female	-5.57589	4.37549	2.17969	-3.77578	-5.59362		
Save	Nanost	tring		_200483200025_R06C01	B cells	Female	-5.91665	4.0189	1.1088	-3.4783	-5.50609		
Save As Text File	Fluidign	n		_200483200025_R07C01	B cells	Female	-5.56902	3.82761	2.02033	-3.29326	-5.34828		
Save As Web Page	Merge Files		_200483200021_R08C01	B cells	Female	-5.65505	4.08722	2.40096	-3.25009	-5.49664			
Close	Plate-b	Plate-based Data		_200483200021_R06C01	B cells	Female	-5.62322	3.91354	2.48638	-4.32497	-5.09504		
Close All	Tagman(g) SAM/BAM file Partek Express Study File BCF/VCF			_200483200025_R01C01	LCLs	Female	-5.11054	4.06107	1.92265	-2.98669	-5.09757		
Close Children				_200483200025_R03C01	LCLs	Female	-5.22815	4.86618	2.30014	-2.97251	-5.33525		
Info				_200483200025_R05C01	LCLs	Female	-4.36911	4.79876	1.23619	-4.53385	-5.13975		
Properties	Zipped	Project.		_200483200025_R02C01	LCLs	Female	-5.35979	4.88132	2.55436	-3.63678	-5.384		
Exit		11	. GSM245211	9_200483200021_R07C01	LCLs	Female	-5.05712	4.72935	2.5945	-3.11181	-5.19711		
		12	. GSM245212	0_200483200021_R05C01	LCLs	Female	-5.01856	5.35611	2.69605	-1.53965	-4.92963		
		Row	s: 12 Column	s: 866839 🖌									
		1.0		(Ť	

Figure 3. Importing the CpG Islands map BED file

• Select the file cpg.bed

The BED file will open as a new spreadsheet.

· Change the spreadsheet name to UCSC CpG Island Annotation and save it

For this region list, you can also calculate the average beta values for the probes in each island per sample and detect differential methylated CpG islands regions. Detailed information on how to get average beta value for each CpG can be found in the *Determining the average values for a region list* section of Starting with a list of genomic regions.

Additional Assistance

If you need additional assistance, please visit our support page to submit a help ticket or find phone numbers for regional support.

