Minimum System Requirements

Partek Genomics Suite (Next Generation Sequencing Studies)

- Windows: 7 SP1 or newer
- Mac: OS 10.11 10.13.4
- Linux: Ubuntu® 12.04, Red Hat® 6, CentOS® 6, or newer
- 64-bit 2GHz quad-core processor
- 16GB of onboard memory
- 500GB of storage available for data and installation
- A graphics card with OpenGL-capable drivers

Partek® Genomics Suite® (Microarray Studies)

- Windows: 7 SP1 or newerMac: OS 10.11 10.13.4
- Linux: Ubuntu[®] 12.04, Red Hat[®] 6, CentOS[®], or newer
- 16GB of onboard memory
- 2GHz processor
- 500MB of storage available for data and installation
- A graphics card with OpenGL-capable drivers

Partek® Pathway™

- Windows: 7 SP1 or newer
- Mac: OS 10.11 10.13.4
- Linux: Ubuntu® 12.04, Red Hat® CentOS®, or newer
- 16GB of onboard memory
- 2GHz processor
- 500MB of storage available for data and installation
- A graphics card with OpenGL-capable drivers

Additional Assistance

If you need additional assistance, please visit our support page to submit a help ticket or find phone numbers for regional support.