Client Computer Connection to License Server

This guide is specific to the client computer connection instructions to a floating concurrent Partek Genomics Suite license server.

Download Partek Genomics Suite

With administrative privileges, download the latest version of Partek Genomics Suite.

Run the Partek Genomics Suite Application

Once the download is completed, start the application by clicking on the **Partek Genomics Suite** icon. The default Partek License Manager window will appear. You will be prompted to provide a license (figure 1).

- 1. Select Add License.
- 2. Select the License server radio button.
- 3. Enter the Server Name and select Add.
- You will need to obtain the server name from your license server administrator.

Partek [®] License Manager Manages your licenses for Partek [®] products 1		Partek*
+ Add License		Partek Licensing Manager - Add License ? ×
Feature	# of Seats Iss	For a stand alone installation, browse to the license file you received from Partek. Otherwise, enter the server details you got from your administrator. Ulicense file Ulicense contents Copy License Content License file directory: C:\Users\rmcdowell\partek_files Browse
Validate Licenses Copy License Information Computer information nmcdowell-pc Host Name nmcdowell-pc Host Ids "b8ac6f9074b4 00ff19fdi Operating System Windows 6.2 x64 (64-bit Copy Information If you currently do not have a license, please contact us) 0441" build)	© License server Server Name: yourservername Port: 3 Add Cancel

Figure 1. Add License

Additional Assistance

If you need additional assistance, please visit our support page to submit a help ticket or find phone numbers for regional support.