## **System Information**

The System information page can be found in Settings>System information (Figure 1).

For new users, Partek Flow provides sample RNA-Seq data available for download. To download the files, select the corresponding data set under Download tutorial data. Depending on the dataset you select, this will automatically create a new project either called Brain vs UHR or iDEA Challenge.

After the data has finished downloading, the samples will automatically be imported into the project. See our RNAseq Tutorial for more information.

Partek Flow details displays the build version currently installed in your server as well as your Root and Library directories, the architecture version and the internet browser you are currently using.

Server details lists additional information about the server that Partek Flow is installed in. Java details lists the version of Java and its configuration as used by Partek Flow. On the other hand, Partek Flow license details lists the licenses available for use in your Partek Flow server.




Figure 1. System information page when viewed by a regular user

## Additional Assistance

If you need additional assistance, please visit our support page to submit a help ticket or find phone numbers for regional support.

