

Post-alignment tools

Post-alignment tools involve tasks that can be performed on aligned data. These are typically used in preparing aligned data for other downstream analyses, such as DNA-Seq or RNA-Seq analysis.

To invoke Post-alignment tools, click on an **Aligned reads** data node (Figure 1). There are three functions available in Post-alignment tools:

- [Filter alignments](#)
- [Convert alignments to unaligned reads](#)
- [Combine alignments](#)
- [Deduplicate UMIs](#)
- [Downscale alignments](#)




Figure 2. Invoking Post-alignment tools

Additional Assistance

If you need additional assistance, please visit [our support page](#) to submit a help ticket or find phone numbers for regional support.



Your Rating:  Results:  43 rates