

Licensing

To access the *Licensing* page, select the avatar in the top right corner, choose **Settings** and select **Licensing** on the left (Figure 1). This page can be accessed by administrative Partek® Flow® users only.

This page displays the status of all license features and toolkits, including their expiration dates. If a license needs to be updated, copy and paste the new license key into the *License key* box and select the **Update** button. For your convenience, the *Partek product* and *Host ID* information is displayed on the screen. If you need a new license key, this information can be copied and sent to licensing@partek.com.

For additional information about the licensing features and toolkits, please see [here](#). Alternatively, please contact your account manager.

Partek Flow

Queue Projects

Administrator

Settings

Help

Contact us

Profile

Log out

Home > Settings > Licensing

Personal settings

My profile

My preferences

System settings

System information

System preferences

LDAP configuration

Partek Flow components

Filter management

Library file management

Option set management

Task management

Pipeline management

Access management

User management

Group management

Licensing

Directory permissions

Access control log

Failed logins

Partek Flow usage

System queue

System resources

Usage report

Current license features

Flow 1000 users, expires 31 December 2017

HPC nodes 1000 nodes, expires 31 December 2017

RNA-seq expires 31 December 2017

User disk quotas expires 31 December 2017

Dedicating licenses to groups expires 31 December 2017

Usage reports expires 31 December 2017

Microarray toolkit expires 31 December 2017

Option sets from command lines expires 31 December 2017

Reprioritizing tasks expires 31 December 2017

REST API expires 31 December 2017

Fluidigm single cell toolkit expires 31 December 2017

User-added tasks expires 31 December 2017

ChIP-seq toolkit expires 31 December 2017

Update license

Copy the license key you received from Partek and paste into the field below then click **Update**.

License key

Update

If you do not have the key, please copy and paste the product and host ID below into an email and send it to licensing@partek.com.

Partek product Flow 6.0

Host ID 0242ac110017,02420a005502

Figure 1. Licensing page

Additional Assistance

If you need additional assistance, please visit [our support page](#) to submit a help ticket or find phone numbers for regional support.



Your Rating:      Results:      36 rates